

Biliary Stents

Single Use

Features:

- Various Types and Sizes Available
- Tapered and Rounded Ends
- Well Balanced Wall Thickness

- Radiopaque Material
- Flaps or Pigtails

Benefits:

- Fits all Anatomic Cases
- Atraumatic Implantation, Easy Placement
- Maximum Inner Lumen with Highest Form Stability and Kink Resistance
- Exact Control of Placement
- Safe and Secure Fixation

Order Information:

Biliary Stent Set, straight with bottom hole

Catalog No.	Stent		Guiding Catheter		Pusher	
	French size	Length between flaps	French size	Length	French size	Length
GBS-01-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm
GBS-01-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm
GBS-01-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm
GBS-01-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm
GBS-01-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm
GBS-01-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm
GBS-01-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm
GBS-01-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm
GBS-01-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm
GBS-01-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm
GBS-01-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm
GBS-01-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm
GBS-01-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm
GBS-01-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm
GBS-01-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm
GBS-01-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm
GBS-01-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm
GBS-01-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm
GBS-01-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm
GBS-01-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm

All „Biliary Stent Sets, straight with bottom hole“, in 7 French include the stent, a pusher and a Teflon® coated guide wire with ø 0.035" and 450 cm length. Stent Sets in 8.5 Fr. and larger include additionally a guiding catheter.

All biliary stents are also available individually. (Catalog No. without "S".)

Manufactured by:

Medi-Globe
Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Cobra Medical BV
Stavangerweg 21-14
9723 JC Groningen, Nederland
Postbus 5153
9700 GD Groningen, Nederland

Contact Nederland

Tel: +31 (0) 50 528 04 98
Fax: +31 (0) 50 528 04 95
Mail: info@cobramedical.nl
www.cobramedical.nl

Contact België

Tel: +32 (0) 70 66 05 70
Fax: +32 (0) 70 66 05 72
Mail: info@cobramedical.be
www.cobramedical.be

Biliary Stent Set, straight with side ports

Catalog No.	Stent		Guiding Catheter		Pusher	
	French size	Length between flaps	French size	Length	French size	Length
GBS-02-07-004S	7.0	4 cm	N/A	N/A	7.0	220 cm
GBS-02-07-006S	7.0	6 cm	N/A	N/A	7.0	220 cm
GBS-02-07-008S	7.0	8 cm	N/A	N/A	7.0	220 cm
GBS-02-07-011S	7.0	11 cm	N/A	N/A	7.0	220 cm
GBS-02-07-016S	7.0	16 cm	N/A	N/A	7.0	220 cm
GBS-02-08-004S	8.5	4 cm	5.0	320 cm	8.5	160 cm
GBS-02-08-006S	8.5	6 cm	5.0	320 cm	8.5	160 cm
GBS-02-08-008S	8.5	8 cm	5.0	320 cm	8.5	160 cm
GBS-02-08-011S	8.5	11 cm	5.0	320 cm	8.5	160 cm
GBS-02-08-016S	8.5	16 cm	5.0	320 cm	8.5	160 cm
GBS-02-10-004S	10.0	4 cm	6.0	320 cm	10.0	160 cm
GBS-02-10-006S	10.0	6 cm	6.0	320 cm	10.0	160 cm
GBS-02-10-008S	10.0	8 cm	6.0	320 cm	10.0	160 cm
GBS-02-10-011S	10.0	11 cm	6.0	320 cm	10.0	160 cm
GBS-02-10-016S	10.0	16 cm	6.0	320 cm	10.0	160 cm
GBS-02-11-004S	11.5	4 cm	6.0	320 cm	11.5	160 cm
GBS-02-11-008S	11.5	6 cm	6.0	320 cm	11.5	160 cm
GBS-02-11-009S	11.5	8 cm	6.0	320 cm	11.5	160 cm
GBS-02-11-011S	11.5	11 cm	6.0	320 cm	11.5	160 cm
GBS-02-11-016S	11.5	16 cm	6.0	320 cm	11.5	160 cm

All "Biliary Stent Sets, straight with side ports", in 7 French include the stent, a pusher and a Teflon® coated guide wire with ø 0.035" and 450 cm length. Stent Sets in 8.5 Fr. and larger include additionally a guiding catheter.

All biliary stents are also available individually. (Catalog No. without "S".)

Biliary Stent Set, Maxi-Flow

Catalog No.	Stent		Guiding Catheter		Pusher	
	French size	Length between flaps	French size	Length	French size	Length
GBS-03-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm
GBS-03-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm
GBS-03-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm
GBS-03-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm
GBS-03-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm
GBS-03-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm
GBS-03-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm
GBS-03-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm
GBS-03-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm
GBS-03-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm
GBS-03-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm
GBS-03-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm
GBS-03-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm
GBS-03-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm
GBS-03-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm
GBS-03-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm
GBS-03-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm
GBS-03-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm
GBS-03-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm
GBS-03-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm

All "Biliary Stent Sets, Maxi Flow", in 7 French include the stent, a pusher and a Teflon® coated guide wire with ø 0.035" and 450 cm length. Stent Sets in 8.5 Fr. and larger include additionally a guiding catheter.

All biliary stents are also available individually. (Catalog No. without "S".)

Biliary Stent Set, Multi-Flap

Catalog No.	Stent		Guiding Catheter		Pusher	
	French size	Length between flaps	French size	Length	French size	Length
GBS-04-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm
GBS-04-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm
GBS-04-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm
GBS-04-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm
GBS-04-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm
GBS-04-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm
GBS-04-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm
GBS-04-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm
GBS-04-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm
GBS-04-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm
GBS-04-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm
GBS-04-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm
GBS-04-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm
GBS-04-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm
GBS-04-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm
GBS-04-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm
GBS-04-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm
GBS-04-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm
GBS-04-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm
GBS-04-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm

All "Biliary Stent Sets, Multi-Flap", in 7 French include the stent, a pusher and a Teflon® coated guide wire with ø 0.035" and 450 cm length. Stent Sets in 8.5 Fr. and larger include additionally a guiding catheter.

All biliary stents are also available individually. (Catalog No. without "S".)

Biliary Stent Set, curved

Catalog No.	Stent		Guiding Catheter		Pusher	
	French size	Length between flaps	French size	Length	French size	Length
GBS-05-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm
GBS-05-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm
GBS-05-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm
GBS-05-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm
GBS-05-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm
GBS-05-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm
GBS-05-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm
GBS-05-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm
GBS-05-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm
GBS-05-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm
GBS-05-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm
GBS-05-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm
GBS-05-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm
GBS-05-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm
GBS-05-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm
GBS-05-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm
GBS-05-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm
GBS-05-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm
GBS-05-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm
GBS-05-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm

All "Biliary Stent Sets, curved", in 7 French include the stent, a pusher and a Teflon® coated guide wire with ø 0.035" and 450 cm length. Stent Sets in 8.5 Fr. and larger include additionally a guiding catheter.

All biliary stents are also available individually. (Catalog No. without "S".)