


Pancreatic Stents

Single Use



Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Cobra Medical BV
Stavangerweg 21-14
9723 JC Groningen, Nederland
Postbus 5153
9700 GD Groningen, Nederland

Contact Nederland

Tel: +31 (0) 50 528 04 98
Fax: +31 (0) 50 528 04 95
Mail: info@cobramedical.nl
www.cobramedical.nl

Contact België

Tel: +32 (0) 70 66 05 70
Fax: +32 (0) 70 66 05 72
Mail: info@cobramedical.be
www.cobramedical.be

Ver. 2.0 - 06/07

 **Medi-Globe**
The Spirit of Care.™


COBRA MEDICAL

Pancreatic Stents

Single Use

Features:

- Various Types and Sizes Available
- Tapered and Rounded Ends
- Radiopaque Material

Benefits:

- > Fits all Anatomic Cases
- > Atraumatic Implantation, Easy Placement
- > Exact Control of Placement

Order Information:

Pancreatic Stent Set, Multi-Flap Pigtail

Catalog No.	Stent		Guide Wire		Pusher	
	French size	Length	Diameter	Length	French size	Length
GBS-08-07-006S	7.0	6 cm	.035" (0.89 mm)	450 cm	7.0	220 cm
GBS-08-07-009S	7.0	9 cm	.035" (0.89 mm)	450 cm	7.0	220 cm
GBS-08-07-012S	7.0	12 cm	.035" (0.89 mm)	450 cm	7.0	220 cm
GBS-08-08-006S	8.5	6 cm	.035" (0.89 mm)	450 cm	8.5	220 cm
GBS-08-08-009S	8.5	9 cm	.035" (0.89 mm)	450 cm	8.5	220 cm
GBS-08-08-012S	8.5	12 cm	.035" (0.89 mm)	450 cm	8.5	220 cm
GBS-08-10-006S	10.0	6 cm	.052" (1.32 mm)	400 cm	10.0	220 cm
GBS-08-10-009S	10.0	9 cm	.052" (1.32 mm)	400 cm	10.0	220 cm
GBS-08-10-012S	10.0	12 cm	.052" (1.32 mm)	400 cm	10.0	220 cm
GBS-08-11-006S	11.5	6 cm	.052" (1.32 mm)	400 cm	11.5	220 cm
GBS-08-11-009S	11.5	9 cm	.052" (1.32 mm)	400 cm	11.5	220 cm
GBS-08-11-012S	11.5	12 cm	.052" (1.32 mm)	400 cm	11.5	220 cm

All „Pancreatic Stent Sets, Multi-Flap Pigtail“, include the stent, a pusher and a Teflon® coated guide wire.
All pancreatic stents are also available individually. (Catalog No. without “S”.)

Pancreatic Stent Set, S-Shaped

Catalog No.	Stent		Guiding Catheter		Pusher	
	French size	Length straightened	French size	Length	French size	Length
GBS-09-05-005S	5.0	5 cm	N/A	N/A	5.0	220 cm
GBS-09-05-007S	5.0	7 cm	N/A	N/A	5.0	220 cm
GBS-09-05-009S	5.0	9 cm	N/A	N/A	5.0	220 cm
GBS-09-05-012S	5.0	12 cm	N/A	N/A	5.0	220 cm
GBS-09-05-015S	5.0	15 cm	N/A	N/A	5.0	220 cm
GBS-09-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm
GBS-09-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm
GBS-09-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm
GBS-09-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm
GBS-09-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm
GBS-09-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm
GBS-09-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm
GBS-09-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm
GBS-09-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm
GBS-09-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm

All “Pancreatic Stent Sets, S-Shape”, include the stent, a pusher and a Teflon® coated guide wire. Stent Sets in 10.0 Fr. include additionally a guiding catheter.

All pancreatic stents are also available individually. (Catalog No. without “S”.)

Pancreatic Stent, straight with radial holes

Catalog-No.	Stent	
	French size	Length
GBS-10-07-005	7.0	5 cm
GBS-10-07-007	7.0	7 cm
GBS-10-07-008	7.0	8 cm
GBS-10-07-009	7.0	9 cm
GBS-10-07-011	7.0	11 cm
GBS-10-07-016	7.0	16 cm

Pancreatic Stent, straight with radial holes

Catalog-No.	Stent	
	French size	Length
GBS-10-08-005	8.5	5 cm
GBS-10-08-007	8.5	7 cm
GBS-10-08-008	8.5	8 cm
GBS-10-08-009	8.5	9 cm
GBS-10-08-011	8.5	11 cm
GBS-10-08-016	8.5	16 cm

Stent-Introducer-Components

Pusher, colored, with metal ring

Catalog No.	French size	Color	Length
GBS-42-05-220	5.0	grey	220 cm
GBS-42-07-220	7.0	yellow	220 cm
GBS-42-08-220	8.5	orange	220 cm
GBS-42-10-220	10.0	red	220 cm
GBS-42-11-220	11.5	blue	220 cm
GBS-42-08-180	8.5	orange	180 cm
GBS-42-10-180	10.0	red	180 cm
GBS-42-11-180	11.5	blue	180 cm
GBS-42-07-160	7.0	yellow	160 cm
GBS-42-08-160	8.5	orange	160 cm
GBS-42-10-160	10.0	red	160 cm
GBS-42-11-160	11.5	blue	160 cm

Guiding Catheter with 1 metal ring

Catalog No.	French size	Length
GBS-41-05-220	5.0	220 cm
GBS-41-05-320	5.0	320 cm
GBS-41-06-220	6.0	220 cm
GBS-41-06-320	6.0	320 cm

Guiding Catheter with 4 metal rings

GBS-43-06-220	6.0	220 cm
GBS-43-06-320	6.0	320 cm

Guide Wire with Teflon® coating, fixed core

Catalog No.	Diameter	Length
GGW-01-35-450	.035" (0.89 mm)	450 cm
GGW-01-52-400	.052" (1.32 mm)	400 cm