


EUS FNA Needle System



EUS FNA Needle System

Endoscopic Ultrasound Needle System for Fine Needle Aspiration

Features:

- Sand-Blasted Needle Tip
- Graduated Needle Adjustor
- Metal Spiral and Aluminum Handle
- Single Use Needle with Stylet

Benefits:

- Optimal Visualization of the Needle
- Variable Needle Projection (0 – 12 cm)
- Spiral and Handle Autoclavable and Reusable
- Safety and Convenience

Order Information:

Catalog No.	Description	Outer Diameter
22 Gauge System		
222280	Aspiration Needle System for Ultrasonic FNA Type "Hancke-Vilmann" consists of Stainless Steel Needle (22 Gauge) with Stylet, Metal Spiral (ø 1.8 mm) and Handle, Sterile Packed, for Pentax Endoscope FG-32UA	1.8 mm
222282	Aspiration Needle for Set 222280, Single Use	
222284	Metal Spiral for Set 222280, Reusable	1.8 mm
22 Gauge System		
222480	Aspiration Needle System for Ultrasonic FNA Type "Hancke-Vilmann" consists of Stainless Steel Needle (22 Gauge) with Stylet, Metal Spiral (ø 1.8 mm) and Handle, Sterile Packed, for Pentax Endoscopes FG-34UX and FG-36UX	1.8 mm
222482	Aspiration Needle for Set 222480, Single Use	
222484	Metal Spiral for Set 222480, Reusable	1.8 mm
19 Gauge System		
223280	Aspiration Needle System for Ultrasonic FNA Type "Hancke-Vilmann" consists of Stainless Steel Needle (19 Gauge) with Stylet, Metal Spiral (ø 2.1 mm) and Handle, Sterile Packed, for Pentax Endoscope FG-38UX	2.1 mm
223282	Aspiration Needle for Set 223280, Single Use	
223284	Metal Spiral for Set 223280, Reusable	2.1 mm
100202	Aluminum Handle for Ultrasound Needle System, for Pentax Endoscopes FG-32UA, FG-34UX, FG-36UX and FG-38UX	

Aluminum Handle



Needle with stylet are for single patient use only. Metal spiral and handle are reusable. Each is sold separately. To ensure compatibility with various marketed Ultrasonic Endoscope types, please contact us.

Manufactured by:

Medi-Globe
Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by: