

# KIMBERLY-CLARK\*

## EASI-LAV\* Closed Gastric Lavage System

### CHAR-FLO\* Activated Charcoal



The KIMBERLY-CLARK\* EASI-LAV\* and CHAR-FLO\* Systems combine to rapidly and effectively cleanse the stomach of patients suffering from accidental poisoning, drug overdoses or gastrointestinal bleeding. Used by both emergency personnel and gastrointestinal lab clinicians, the EASI-LAV and CHAR-FLO systems are available in sizes from pediatric to adult, and can complete gastric lavage and charcoal delivery quickly, reducing patient's exposure to toxins.

In addition to rapid and efficient patient treatment, the KIMBERLY-CLARK EASI-LAV and CHAR-FLO active closed systems help prevent caregivers from exposure to patient's body fluids. Additional features include:

- **Dual plungers doubling the fluid exchange volume to reduce lavage time and ensure effective particulate removal**
- **CHAR-FLO attaches directly to a plunger, helps eliminate risk of accidental disconnections that cause "blackouts"**
- **Complete kits or single tubes**

From medical devices to surgical apparel – everything we do, everything we make, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

KIMBERLY-CLARK\*  
EASI-LAV\* Closed Gastric Lavage System  
CHAR-FLO\* Activated Charcoal



*Dual plungers doubling the fluid exchange volume to reduce lavage time and ensure effective particulate removal*



*CHAR-FLO attaches directly to a plunger, helps eliminate risk of accidental disconnections that cause "blackouts"*

**KIMBERLY-CLARK\* EASI-LAV\* Closed Gastric Lavage System**

Stock #	Size	Description	Packaging
1518	18 Fr	EASI-LAV with Gastric Tube (with vent)	5/Case
1524	24 Fr	EASI-LAV with Gastric Tube	5/Case
1528	28 Fr	EASI-LAV with Gastric Tube	5/Case
1534	34 Fr	EASI-LAV with Gastric Tube	5/Case
15341	34 Fr	EASI-LAV with Gastric Tube	5/Case
1540	40 Fr	EASI-LAV with Gastric Tube (with 1503 CHAR-FLO*)	5/Case
1501	N/A	Waste Bags	10/Case

**KIMBERLY-CLARK\* EASI-TUBE Gastric Tube**

Stock #	Size	Description	Packaging
15018	18 Fr	Gastric Tube (with vent)	5/Case
15024	24 Fr	Gastric Tube	5/Case
15028	28 Fr	Gastric Tube	5/Case
15034	34 Fr	Gastric Tube	5/Case
15040	40 Fr	Gastric Tube	5/Case

**KIMBERLY-CLARK\* CHAR-FLO\* Activated Charcoal**

Stock #	Size	Description	Packaging
1502	25g/120ml	Charcoal with Sorbitol	12/Case
1504	25g/120ml	Charcoal with Aqueous Base	12/Case
1503	50g/240ml	Charcoal with Sorbitol	12/Case
1505	50g/240ml	Charcoal with Aqueous Base	12/Case



**Commitment to Excellence**

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-528-5591 in the United States, +801 572 6800 internationally, or visit our web site at [www.kchealthcare.com](http://www.kchealthcare.com).

\* Registered Trademark of Kimberly-Clark Corp. Roswell, GA 30076 or its affiliates. Ballard Medical Products, a wholly-owned subsidiary of Kimberly-Clark Corporation, Draper, Utah 84020 USA. © 2001 KCC. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Kimberly-Clark's Neenah Paper.

Z0381 U3324

Protection. For life.

