

E.R.C.P. Catheters

Single Use

E.R.C.P. Catheters


- "Pre-Curved" Catheter with Memory Effect
Easier Positioning to the Papilla and Possible Manual Tip Manipulation for Individual Anatomy
- Guide Wire Compatibility
Contrast Injection while Guide Wire is in Place
- Metal Tips / Metal Rings
Maximum X-Ray Visualization
- Color Graduation
Macroscopically Controlled Cannulation Depth



E.R.C.P. Catheters

Single Use

E.R.C.P. Catheters



Compatible to .035“ Guide Wire

E.R.C.P. Catheter – “Pre-Curved” with Color Graduation and Metal Tips / Metal Rings, Removable Nitinol Stylet, sterile packed

| Catalog No. | Tip Shape | Catheter Outer Diameter |
|---------------|-----------------------|-------------------------|
| GER-01-18-035 | Standard Tip | 1.8 mm |
| GER-02-18-035 | 5 mm Tapered Tip | 1.8 mm |
| GER-04-18-035 | Metal Cone Tip | 1.8 mm |
| GER-05-18-035 | Round Metal Tip | 1.8 mm |
| GER-06-18-035 | Bullet Type Metal Tip | 1.8 mm |
| GER-07-18-035 | Metal Ball Tip | 1.8 mm |

Compatible to .021“ Guide Wire

E.R.C.P. Catheter – “Pre-Curved” with Color Graduation and Metal Tips / Metal Rings, Removable Nitinol Stylet, sterile packed

| Catalog No. | Tip Shape | Catheter Outer Diameter |
|---------------|-------------------|-------------------------|
| GER-03-18-021 | 10 mm Tapered Tip | 1.8 mm |

Distributed by: