



Injectra

Injection Needles – Single Use

Injectra

- Top Down Injection
User-friendly Injection of Medicine
- Removable Safety Clip
Prevents Needle Projection during Accessory Channel Advancement
- Automatic Needle Retraction by Pressing a Release Button
Secures the Needle in the Plastic Sheath during Instrument Removal
- Needle Lock Mechanism
Stable & Controlled Puncture in Targeted Tissue and Vessels



 **Medi-Globe**[®]
The Spirit of Care.[™]

Bronchoscopy

Injectra – Injection needle, sterile packed

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length
GSN-04-23-120	19 Gauge / 1.0 mm	6.0 mm	2.3 mm	120 cm
GSN-03-19-120	22 Gauge / 0.7 mm	6.0 mm	1.9 mm	120 cm

Gastroscopy

Injectra – Injection needle, sterile packed

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length
GSN-02-23-180	19 Gauge / 1.0 mm	4.0 mm	2.3 mm	180 cm
GSN-04-23-160	19 Gauge / 1.0 mm	6.0 mm	2.3 mm	160 cm
GSN-04-23-180	19 Gauge / 1.0 mm	6.0 mm	2.3 mm	180 cm
GSN-01-19-180	22 Gauge / 0.7 mm	4.0 mm	1.9 mm	180 cm
GSN-03-19-160	22 Gauge / 0.7 mm	6.0 mm	1.9 mm	160 cm
GSN-03-19-180	22 Gauge / 0.7 mm	6.0 mm	1.9 mm	180 cm
GSN-05-19-180	25 Gauge / 0.5 mm	4.0 mm	1.9 mm	180 cm
GSN-06-19-180	25 Gauge / 0.5 mm	6.0 mm	1.9 mm	180 cm

Colonoscopy

Injectra – Injection needle, sterile packed

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length
GSN-02-23-230	19 Gauge / 1.0 mm	4.0 mm	2.3 mm	230 cm
GSN-04-23-230	19 Gauge / 1.0 mm	6.0 mm	2.3 mm	230 cm
GSN-01-19-230	22 Gauge / 0.7 mm	4.0 mm	1.9 mm	230 cm
GSN-03-19-230	22 Gauge / 0.7 mm	6.0 mm	1.9 mm	230 cm
GSN-05-19-230	25 Gauge / 0.5 mm	4.0 mm	1.9 mm	230 cm
GSN-06-19-230	25 Gauge / 0.5 mm	6.0 mm	1.9 mm	230 cm