

Cytology Brushes

Single Use

Cytology Brushes

- Brush with Ellipsoid Tip
Efficient Tube Closure to Secure Smear Yield
- Set Versions including Brush and Scissor
Safe and User-friendly Specimen Handling
- Sophexa Easy Brush
Guide Wire Compatible Cytology Brush for Smear Acquisition from the Bile Duct



Cytology Brushes

Single Use

Bronchoscopy

Cytology Brush – with Ellipsoid Tip and Ergonomic Handle

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-01-18-120	Cytology Brush	120 cm	1.8 mm	2.0 mm
GCB-01-24-120	Cytology Brush	120 cm	2.45 mm	2.8 mm
GCB-02-18-120	Set incl. Cytology Brush & Scissor	120 cm	1.8 mm	2.0 mm
GCB-02-24-120	Set incl. Cytology Brush & Scissor	120 cm	2.45 mm	2.8 mm

Gastroscopy

Cytology Brush – with Ellipsoid Tip and Ergonomic Handle

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-01-18-180	Cytology Brush	180 cm	1.8 mm	2.0 mm
GCB-01-24-180	Cytology Brush	180 cm	2.45 mm	2.8 mm
GCB-02-18-180	Set incl. Cytology Brush & Scissor	180 cm	1.8 mm	2.0 mm
GCB-02-24-180	Set incl. Cytology Brush & Scissor	180 cm	2.45 mm	2.8 mm

Colonoscopy

Cytology Brush – with Ellipsoid Tip and Ergonomic Handle

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-01-18-230	Cytology Brush	230 cm	1.8 mm	2.0 mm
GCB-01-24-230	Cytology Brush	230 cm	2.45 mm	2.8 mm
GCB-02-18-230	Set incl. Cytology Brush & Scissor	230 cm	1.8 mm	2.0 mm
GCB-02-24-230	Set incl. Cytology Brush & Scissor	230 cm	2.45 mm	2.8 mm

Sophexa Easy Brush

Guide Wire Compatible Cytology Brush for Smear Acquisition mainly from the Bile Duct.

Cytology Brush – with Ellipsoid Tip and Ergonomic Handle

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-03-24-230	Cytology Brush	230 cm	2.45 mm	2.8 mm