


Pancreatic Stents

Single Use



Distributed by:

Ver. 3.0 (2010-10)

Manufactured by:

Pancreatic Stents

Single Use

Pancreatic Stents

- Well Balanced Wall Diameter
Maximum Inner Lumen with Highest Form Stability & Kink Resistance
- Tapered & Round Ends
Atraumatic Implantation
- Radiopaque Material
Exact Control of Placement



Pancreatic Stents


Medi-Globe offers the pancreatic stents as single component (catalog no. without "S") and additionally as integrated part of a set. Depending on the stent diameter and shape sets include the following components:

7 French Stents and All Pigtail / Double Pigtail Stents

- Stent
- Pusher
- Teflon Coated Stainless Steel Guide Wire


Straight Stents in Diameter Larger than 7 French

- Stent
- Pusher
- Teflon Coated Stainless Steel Guide Wire
- Guiding Catheter



Pancreatic Stent Set Multi-Flap with Pigtail

| Catalog No. | Stent | | Pusher | | Guide Wire | |
|----------------|-------------|------------------------|-------------|--------|-----------------|--------|
| | French Size | Length between Pigtail | French Size | Length | Diameter | Length |
| GBS-08-07-006S | 7.0 | 6 cm | 7.0 | 220 cm | .035" (0.89 mm) | 450 cm |
| GBS-08-07-009S | 7.0 | 9 cm | 7.0 | 220 cm | .035" (0.89 mm) | 450 cm |
| GBS-08-07-012S | 7.0 | 12 cm | 7.0 | 220 cm | .035" (0.89 mm) | 450 cm |
| GBS-08-08-006S | 8.5 | 6 cm | 8.5 | 220 cm | .035" (0.89 mm) | 450 cm |
| GBS-08-08-009S | 8.5 | 9 cm | 8.5 | 220 cm | .035" (0.89 mm) | 450 cm |
| GBS-08-08-012S | 8.5 | 12 cm | 8.5 | 220 cm | .035" (0.89 mm) | 450 cm |
| GBS-08-10-006S | 10.0 | 6 cm | 10.0 | 220 cm | .052" (1.32 mm) | 400 cm |
| GBS-08-10-009S | 10.0 | 9 cm | 10.0 | 220 cm | .052" (1.32 mm) | 400 cm |
| GBS-08-10-012S | 10.0 | 12 cm | 10.0 | 220 cm | .052" (1.32 mm) | 400 cm |
| GBS-08-11-006S | 11.5 | 6 cm | 11.5 | 220 cm | .052" (1.32 mm) | 400 cm |
| GBS-08-11-009S | 11.5 | 9 cm | 11.5 | 220 cm | .052" (1.32 mm) | 400 cm |
| GBS-08-11-012S | 11.5 | 12 cm | 11.5 | 220 cm | .052" (1.32 mm) | 400 cm |




Pancreatic Stent Straight with Radial Holes

| Catalog No. | French Size | Length between Flap |
|---------------|-------------|---------------------|
| GBS-10-07-005 | 7.0 | 5 cm |
| GBS-10-07-007 | 7.0 | 7 cm |
| GBS-10-07-008 | 7.0 | 8 cm |
| GBS-10-07-009 | 7.0 | 9 cm |
| GBS-10-07-011 | 7.0 | 11 cm |
| GBS-10-07-016 | 7.0 | 16 cm |

Pancreatic Stent Straight with Radial Holes

| Catalog No. | French size | Length between Flap |
|---------------|-------------|---------------------|
| GBS-10-08-005 | 8.5 | 5 cm |
| GBS-10-08-007 | 8.5 | 7 cm |
| GBS-10-08-008 | 8.5 | 8 cm |
| GBS-10-08-009 | 8.5 | 9 cm |
| GBS-10-08-011 | 8.5 | 11 cm |
| GBS-10-08-016 | 8.5 | 16 cm |



Pancreatic Stent Set S - Shape

| Catalog No. | Stent | | Guiding Catheter | | Pusher | | Guide Wire | |
|----------------|-------------|--------|------------------|--------|-------------|--------|-----------------|--------|
| | French Size | Length | French Size | Length | French Size | Length | Diameter | Length |
| GBS-09-05-005S | 5.0 | 5 cm | N/A | N/A | 5.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-05-007S | 5.0 | 7 cm | N/A | N/A | 5.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-05-009S | 5.0 | 9 cm | N/A | N/A | 5.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-05-012S | 5.0 | 12 cm | N/A | N/A | 5.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-05-015S | 5.0 | 15 cm | N/A | N/A | 5.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-07-005S | 7.0 | 5 cm | N/A | N/A | 7.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-07-007S | 7.0 | 7 cm | N/A | N/A | 7.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-07-009S | 7.0 | 9 cm | N/A | N/A | 7.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-07-012S | 7.0 | 12 cm | N/A | N/A | 7.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-07-015S | 7.0 | 15 cm | N/A | N/A | 7.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-10-005S | 10.0 | 5 cm | 6.0 | 320 cm | 10.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-10-007S | 10.0 | 7 cm | 6.0 | 320 cm | 10.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-10-009S | 10.0 | 9 cm | 6.0 | 320 cm | 10.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-10-012S | 10.0 | 12 cm | 6.0 | 320 cm | 10.0 | 220 cm | .035" (.089 mm) | 450 cm |
| GBS-09-10-015S | 10.0 | 15 cm | 6.0 | 320 cm | 10.0 | 220 cm | .035" (.089 mm) | 450 cm |

Guiding Catheters

Guiding Catheters are available with one and with four metal rings. The metal markers on the 4 ring version are in a distance of 5 cm and allow a pre-measurement of the stricture length and ease the selection of the correct stent.



Guiding Catheter with 1 Metal Ring

| Catalog No. | French Size | Length |
|---------------|-------------|--------|
| GBS-41-05-220 | 5.0 | 220 cm |
| GBS-41-05-320 | 5.0 | 320 cm |
| GBS-41-06-220 | 6.0 | 220 cm |
| GBS-41-06-320 | 6.0 | 320 cm |



Guiding Catheter with 4 Metal Rings

| Catalog No. | French Size | Length |
|---------------|-------------|--------|
| GBS-43-06-220 | 6.0 | 220 cm |
| GBS-43-06-320 | 6.0 | 320 cm |

Pushers

Pushers follow for each diameter a defined color concept to differentiate the related French sizes.



Pusher with 1 Metal Ring

| Catalog No. | French Size | Color Concept | Length |
|---------------|-------------|---------------|--------|
| GBS-42-05-220 | 5.0 | grey | 220 cm |
| GBS-42-07-160 | 7.0 | yellow | 160 cm |
| GBS-42-07-220 | 7.0 | yellow | 220 cm |
| GBS-42-08-160 | 8.5 | orange | 160 cm |
| GBS-42-08-180 | 8.5 | orange | 180 cm |
| GBS-42-08-220 | 8.5 | orange | 220 cm |
| GBS-42-10-160 | 10.0 | red | 160 cm |
| GBS-42-10-180 | 10.0 | red | 180 cm |
| GBS-42-10-220 | 10.0 | red | 220 cm |
| GBS-42-11-160 | 11.5 | blue | 160 cm |
| GBS-42-11-180 | 11.5 | blue | 180 cm |
| GBS-42-11-220 | 11.5 | blue | 220 cm |