


# RotaCut<sup>®</sup>

Sphincterotomes - Single Use



Manufactured by:

**Medi-Globe**

Medi-Globe GmbH  
Medi-Globe-Strasse 1-5  
83101 Achenmühle, Germany  
Phone: ++49 8032 973 379  
Fax: ++49 8032 973 339  
e-mail: sales@medi-globe.de  
www.medi-globe.de

Distributed by:

**Cobra Medical BV**  
Stavangerweg 21-14  
9723 JC Groningen, Nederland  
Postbus 5153  
9700 GD Groningen, Nederland

**Contact Nederland**

Tel: +31 (0) 50 528 04 98  
Fax: +31 (0) 50 528 04 95  
Mail: info@cobramedical.nl  
www.cobramedical.nl

**Contact België**

Tel: +32 (0) 70 66 05 70  
Fax: +32 (0) 70 66 05 72  
Mail: info@cobramedical.be  
www.cobramedical.be

Ver. 2.0 - 06/07

 **Medi-Globe**  
The Spirit of Care.™

 **COBRA MEDICAL**

# RotaCut®

Sphincterotomes - Single Use

## Features:

- 2 & 3 Lumen Versions
- Rotatable Cutting Wire
- Color Markings with Black Rings
- "Pre-curved" Versions
- Ergonomically Shaped Handle

## Benefits:


- Contrast Medium Injection while Guide Wire is in Place
- Precise Positioning
- Macroscopically Controlled Cannulation Depth
- Easier Positioning to the Papilla
- User-Friendly Handling

## Order Information:

### 3-Lumen Guide Wire Sphincterotomes, rotatable:

Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-31-20-020	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-31-20-025	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	25 mm	.035" (0.89 mm)
GSP-31-20-030	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)

### Tips:



### 2-Lumen Guide Wire Sphincterotomes, rotatable:

Catalog Nr.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-21-17-020	Tapered Tip, Tip Length 5 mm	1.7 mm	20 mm	.021" (0.53 mm)
GSP-21-17-030	Tapered Tip, Tip Length 5 mm	1.7 mm	30 mm	.021" (0.53 mm)
GSP-22-17-020	Precut	1.7 mm	20 mm	.021" (0.53 mm)
GSP-22-17-030	Precut	1.7 mm	30 mm	.021" (0.53 mm)
GSP-21-20-020	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-21-20-030	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-24-20-020	Tapered Tip, Tip Length 10 mm	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-24-20-030	Tapered Tip, Tip Length 10 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-23-20-030	Tapered Tip, Tip Length 20 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-22-20-020	Precut	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-22-20-030	Precut	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-34-20-020	B II Sphincterotome	2.55/2.05 mm	Variable	.035" (0.89 mm)

### Standard-Sphincterotomes, 1-Lumen, rotatable:

Catalog No.	Description	Outer Diameter	Cutting Wire Length	Tip Length
GSP-01-17-020	Standard Tip	1.7 mm	20 mm	5.0 mm
GSP-01-17-030	Standard Tip	1.7 mm	30 mm	5.0 mm
GSP-02-17-020	Precut	1.7 mm	20 mm	-
GSP-02-17-030	Precut	1.7 mm	30 mm	-
GSP-03-17-020	Precut with Short Tip	1.7 mm	20 mm	1.0 mm
GSP-03-17-030	Precut with Short Tip	1.7 mm	30 mm	1.0 mm
GSP-04-17-020	Tapered Tip	1.7 mm	20 mm	5.0 mm
GSP-04-17-030	Tapered Tip	1.7 mm	30 mm	5.0 mm
GSP-14-18-030	B II Sphincterotome	1.8 mm	Variable	-

### Needle Knife Sphincterotomes, 1-Lumen, non-rotatable:

Catalog No.	Description	Outer Diameter
GSP-11-17-005	Needle Knife Sphincterotome, Needle Projection 0 - 5 mm	1.7 mm
GSP-11-17-010	Needle Knife Sphincterotome, Needle Projection 0 - 10 mm	1.7 mm
GSP-11-17-015	Needle Knife Sphincterotome, Needle Projection 0 - 15 mm	1.7 mm

### Special-Sphincterotomes, non-rotatable:

Catalog No.	Description	Outer Diameter
GSP-12-18-004	Fistolotome; cutting spiral 4 mm; color markings	1.8 mm
GSP-13-17-015	Shark Fin Sphincterotome; cutting wire 15 mm; color markings	1.7 mm

Medi-Globe handles are directly compatible with the Olympus and Microvasive electro-surgical generators. To ensure compatibility with other types of electro-surgical generators, please contact Medi-Globe GmbH or your distributor. We recommend the use of compatible Medi-Globe active cords. You will find the corresponding Medi-Globe active cords in the leaflet "Accessories".