

# Accessories



Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH  
Medi-Globe-Strasse 1-5  
83101 Achenmühle, Germany  
Phone: ++49 8032 973 379  
Fax: ++49 8032 973 339  
e-mail: sales@medi-globe.de  
www.medi-globe.de

Distributed by:

Rev. 01 - 04/06


 **Medi-Globe**<sup>®</sup>  
*The Spirit of Care*<sup>™</sup>

# Spray Catheters

Reusable

## Description:

Coaxial Spray Catheters for staining and anesthetizing, autoclavable



## Order Information:


| Catalog No. | Description             | Spray Head Diameter | Catheter Diameter | Length |
|-------------|-------------------------|---------------------|-------------------|--------|
| 540360      | with Coaxial Spray Head | 2.3 mm              | 1.8 mm            | 180 cm |
| 540460      | with Coaxial Spray Head | 2.3 mm              | 1.8 mm            | 230 cm |

# Coagulation Probes

Reusable

## Description:

Monopolar Coagulation Probes with atraumatic dome tip, autoclavable



## Order Information:

| Catalog No. | Diameter | Length |
|-------------|----------|--------|
| 730360      | 1.8 mm   | 180 cm |
| 731360      | 2.4 mm   | 180 cm |
| 731460      | 2.4 mm   | 230 cm |

# Active Cords

Reusable

## Description:

Active Cords for connecting Medi-Globe products to generators of various manufacturers

## Order Information:



| Catalog No.   | Description                              | Length |
|---------------|--|--------|
| GAC-51-04-300 | Active Cord for Erbe 4 mm                | 300 cm |
| GAC-51-05-300 | Active Cord for Erbe 5 mm                | 300 cm |
| GAC-52-04-300 | Active Cord for Martin / Berchthold 4 mm | 250 cm |
| GAC-53-08-300 | Active Cord for Valleylab 8 mm           | 300 cm |
| GAC-54-06-300 | Active Cord for Olympus 6 mm             | 300 cm |

# Pill Splitter, Pill Crusher & Combination


Reusable

## Features:

- Dishwasher-Safe
- Storage Container
- Surgical Stainless Steel Cutting Blade
- Unique, Ergonomic Design

## Benefits:


- Easy Cleaning
- Safe Storage of Pills
- Precise Application and Durability
- Convenience and Stability



Pill Splitter



Pill Crusher



Pill Splitter / Pill Crusher-Combination

## Order Information:

| Catalog No. | Description                           |
|-------------|---------------------------------------|
| 10401       | Pill Splitter                         |
| 10402       | Pill Crusher                          |
| 10403       | Pill Splitter / Crusher - Combination |