

Bavarian Wires

Guide Wires



 **Medi-Globe**[®]
The Spirit of Care.[™]

 **COBRA MEDICAL**

Bavarian Wires

Guide Wires

Features:

- Nitinol Versions
- Hydrophilic Coated Versions
- Teflon® Coated
- Two Color Markings and Radiopaque Tip
- Single Use or Autoclavable

Benefits:

- Kink Resistant
- Easy, Low Friction Manipulation
- Usable during Sphincterotomy
- Enhanced Visualization
- Safety and Convenience or Price Reduction by Multi-Use

Order Information:

Catalog No. Single Use	Description Teflon® coated, stainless steel	Diameter	Length
---------------------------	--	----------	--------

GGW-01-21-400	Fixed Core	.021" (0.53 mm)	400 cm
GGW-01-21-450	Fixed Core	.021" (0.53 mm)	450 cm
GGW-01-35-400	Fixed Core	.035" (0.89 mm)	400 cm
GGW-01-35-450	Fixed Core	.035" (0.89 mm)	450 cm
GGW-01-47-400	Fixed Core	.047" (1.20 mm)	400 cm
GGW-01-47-450	Fixed Core	.047" (1.20 mm)	450 cm
GGW-01-52-400	Fixed Core	.052" (1.32 mm)	400 cm
GGW-01-52-450	Fixed Core	.052" (1.32 mm)	450 cm
GGW-02-35-400	Movable Core	.035" (0.89 mm)	400 cm
GGW-02-35-450	Movable Core	.035" (0.89 mm)	450 cm

Catalog No. Autoclavable	Description Bavarian Wires – Hydrophilic Teflon® coated Nitinol Wire, two-color markings	Diameter	Length
-----------------------------	---	----------	--------

GGW-03-21-400	"Soft" Tip, 50 mm tapered radiopaque tip	.021" (0.53 mm)	400 cm
GGW-03-21-450	"Soft" Tip, 50 mm tapered radiopaque tip	.021" (0.53 mm)	450 cm
GGW-03-25-450	"Soft" Tip, 50 mm tapered radiopaque tip	.025" (0.64 mm)	450 cm
GGW-03-35-400	"Soft" Tip, 50 mm tapered radiopaque tip	.035" (0.89 mm)	400 cm
GGW-03-35-450	"Soft" Tip, 50 mm tapered radiopaque tip	.035" (0.89 mm)	450 cm
GGW-04-35-400	"Standard" Tip, 50 mm tapered radiopaque tip	.035" (0.89 mm)	400 cm
GGW-04-35-450	"Standard" Tip, 50 mm tapered radiopaque tip	.035" (0.89 mm)	450 cm

Catalog No. Autoclavable	Description Bavarian Wires – for Short Wire Systems, Hydrophilic Teflon® coated Nitinol Wire, two-color markings	Diameter	Length
-----------------------------	---	----------	--------

GGW-03-21-260	"Soft" Tip, 50 mm tapered radiopaque tip	.021" (0.53 mm)	260 cm
GGW-03-25-260	"Soft" Tip, 50 mm tapered radiopaque tip	.025" (0.64 mm)	260 cm
GGW-03-35-260	"Soft" Tip, 50 mm tapered radiopaque tip	.035" (0.89 mm)	260 cm

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Cobra Medical BV
Stavangerweg 21-14
9723 JC Groningen, Nederland
Postbus 5153
9700 GD Groningen, Nederland

Contact Nederland

Tel: +31 (0) 50 528 04 98
Fax: +31 (0) 50 528 04 95
Mail: info@cobramedical.nl
www.cobramedical.nl

Contact België

Tel: +32 (0) 70 66 05 70
Fax: +32 (0) 70 66 05 72
Mail: info@cobramedical.be
www.cobramedical.be