


Nasal Biliary Drainages

Single Use



Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Cobra Medical BV
Stavangerweg 21-14
9723 JC Groningen, Nederland
Postbus 5153
9700 GD Groningen, Nederland

Contact Nederland

Tel: +31 (0) 50 528 04 98
Fax: +31 (0) 50 528 04 95
Mail: info@cobramedical.nl
www.cobramedical.nl

Contact België

Tel: +32 (0) 70 66 05 70
Fax: +32 (0) 70 66 05 72
Mail: info@cobramedical.be
www.cobramedical.be

Ver. 2.0 - 06/07

 **Medi-Globe**
The Spirit of Care.™

 **COBRA MEDICAL**

Nasal Biliary Drainages

Single Use


Features:

- Various Types and Sizes Available
- Tapered and Rounded Ends
- Well-Balanced Wall Thickness
- Radiopaque Material
- Single Use

Benefits:

- Fits all Anatomic Cases
- Atraumatic Implantation, Easy Placement
- Maximum Inner Lumen with Guaranteed Form Stability and Kink Resistance
- Exact Positioning
- Safety and Convenience

Nasal Biliary Drainage Sets:



Order Information:

Nasal Biliary Drainage Sets

Catalog No.	Description	Drainage Catheter		Nasal Tube Length	Drainage Connecting Tube Length	Recommended Guide Wire	
		French size	Length			Diameter	Length
GND-01-05-290S	Pigtail	5.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-01-07-290S	Pigtail	7.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-01-10-290S	Pigtail	10.0	290 cm	50 cm	135 cm	.052" (1.32 mm)	400 cm
GND-02-05-290S	Standard	5.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-02-07-290S	Standard	7.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-02-10-290S	Standard	10.0	290 cm	50 cm	135 cm	.052" (1.32 mm)	400 cm
GND-03-05-290S	Standard with Pigtail	5.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-03-07-290S	Standard with Pigtail	7.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-03-10-290S	Standard with Pigtail	10.0	290 cm	50 cm	135 cm	.052" (1.32 mm)	400 cm

All Nasal Biliary Drainage Sets include a Drainage Catheter, Nasal Tube, Drainage Connecting Tube, Injection Connector and a Teflon® coated Guide Wire.

All Nasal Biliary Drainage Catheters are also available individually (catalog no. without „S“).

Nasal Tube

Catalog No.	French size	Length
GND-41-05-050	5.0	50 cm
GND-41-07-050	7.0	50 cm
GND-41-10-050	10.0	50 cm




Injection Connector

Catalog No.	French size
GND-43-05-000	5.0
GND-43-07-000	7.0
GND-43-10-000	10.0

Guide Wire with Teflon®-coating, fixed core

Catalog No.	Diameter	Length
GGW-01-35-400	.035" (0.89 mm)	400 cm
GGW-01-52-400	.052" (1.32 mm)	400 cm



Drainage Connecting Tube

Catalog No.	Length
GND-42-03-135	135 cm