

StoneBuster®

Flushable Mechanical Lithotripsy System

FLUSHABLE METAL SPIRAL

- › Cannulation via the New Flushable Metal Spiral and Basket - No Exchanges or Additional Components Required!
- › Direct Introduction of Contrast Medium
- › System Components Available in Reusable or Single Use Versions

STONEBUSTER® HANDLE

- › Precise Manipulation of the Basket
- › Unique Twist and Pull Action Delivers High-Basket Tension for Difficult Biliary Stones
- › User Friendly - Quick & Easy to Assemble and Clean

 **Medi-Globe®**
The Spirit of Care.™

 **COBRA MEDICAL**

StoneBuster®

Flushable Mechanical Lithotripsy System

Lithotripsy System Set - Flushable

GML-99-02-000S Mechanical Lithotripsy Set with Flushable **Single Use** Baskets including a Mechanical Lithotripsy Handle, two Flushable Metal Spirals, one Emergency Spiral, one Instrument Holder (100601) and three Single Use Baskets (GML-02-26-425, GML-02-26-430 and GML-02-26-630)

GML-99-52-000S Mechanical Lithotripsy Set with Flushable **Reusable** Baskets Including a Mechanical Lithotripsy Handle, two Flushable Metal Spirals, one Emergency Spiral, one Instrument Holder (100601) and two Reusable Baskets (GML-52-26-425 and GML-52-26-430)

Catalog No.	Description	Outer ø	Basket ø	Minimum Working Channel ø
-------------	-------------	---------	----------	---------------------------

Lithotripsy Basket Sets, pre-assembled-Single Use Basket with Flushable Metal Spiral and Instrument Holder, sterile packed

GML-11-26-420	4-Wire Basket	2.6 mm	20 mm	2.8 mm
GML-11-26-425	4-Wire Basket	2.6 mm	25 mm	2.8 mm
GML-11-26-430	4-Wire Basket	2.6 mm	30 mm	2.8 mm
GML-11-26-620	6-Wire Basket	2.6 mm	20 mm	2.8 mm
GML-11-26-625	6-Wire Basket	2.6 mm	25 mm	2.8 mm
GML-11-26-630	6-Wire Basket	2.6 mm	30 mm	2.8 mm

Single Use Lithotripsy Basket (only), Length 280 cm, sterile packed

GML-02-26-420	4-Wire Basket	2.6 mm	20 mm	2.8 mm
GML-02-26-425	4-Wire Basket	2.6 mm	25 mm	2.8 mm
GML-02-26-430	4-Wire Basket	2.6 mm	30 mm	2.8 mm
GML-02-26-620	6-Wire Basket	2.6 mm	20 mm	2.8 mm
GML-02-26-625	6-Wire Basket	2.6 mm	25 mm	2.8 mm
GML-02-26-630	6-Wire Basket	2.6 mm	30 mm	2.8 mm

Reusable Lithotripsy Basket (only), Length 280 cm

GML-52-26-420	4-Wire Basket	2.6 mm	20 mm	2.8 mm
GML-52-26-425	4-Wire Basket	2.6 mm	25 mm	2.8 mm
GML-52-26-430	4-Wire Basket	2.6 mm	30 mm	2.8 mm
GML-52-26-620	6-Wire Basket	2.6 mm	20 mm	2.8 mm
GML-52-26-625	6-Wire Basket	2.6 mm	25 mm	2.8 mm
GML-52-26-630	6-Wire Basket	2.6 mm	30 mm	2.8 mm

Lithotripsy System Components

Catalog No.	Description	Outer ø	Minimum Working Channel ø
GML-99-01-000	Patented Lithotriptor Handle		
GML-91-26-180	Flushable Metal Spiral, Length 180 cm	2.6 mm	2.8 mm
GML-90-34-080	Emergency Spiral, Length 80 cm	3.4 mm	-
100601	Instrument Holder (for Flushable Baskets)		

Manufactured by:



Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Cobra Medical BV
Stavangerweg 21-14
9723 JC Groningen, Nederland
Postbus 5153
9700 GD Groningen, Nederland

Contact Nederland

Tel: +31 (0) 50 528 04 98
Fax: +31 (0) 50 528 04 95
Mail: info@cobramedical.nl
www.cobramedical.nl

Contact België

Tel: +32 (0) 70 66 05 70
Fax: +32 (0) 70 66 05 72
Mail: info@cobramedical.be
www.cobramedical.be