




Stone Extraction Balloons

Single Use



2 & 3-Lumen
Versions



 **Medi-Globe**[®]
The Spirit of Care.™

 **COBRA MEDICAL**

Stone Extraction Balloons

Single Use

Features:

- 2 & 3 Lumen Versions
- Two-Step Balloon Inflation
- Metal Markings at the Distal and Proximal End of the Balloon
- One Centimeter Markings (30 cm)
- Removable Stylet (1.75 m)

Benefits:

- Guide Wire and Contrast Medium Compatibility
- Two Defined Balloon Diameters per Size
- Exact Positioning
- Macroscopically Controlled Penetration Depth
- Kink Free Advancement through the Working Channel

Order Information:

Catalog No.	French size	Length	∅ Balloon	Max. Air Inflation	Max. Liquid Inflation	Accepts Guide Wire
3-Lumen Stone Extraction Balloon, single use						
GEB-02-25-035	7.5	200 cm	14-16 mm	4.5 ml	2.5 ml	.035" (0.89 mm)
2-Lumen Stone Extraction Balloons, single use						
GEB-01-17-025	5.0	200 cm	9-12 mm	3.0 ml	1.5 ml	.021" (0.53 mm)
GEB-01-23-035	7.0	200 cm	14-16 mm	4.5 ml	2.5 ml	.035" (0.89 mm)

Caution: Balloon Material contains Latex.

Manufactured by:

Medi-Globe
 Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by:

Cobra Medical BV
 Stavangerweg 21-14
 9723 JC Groningen, Nederland
 Postbus 5153
 9700 GD Groningen, Nederland

Contact Nederland

Tel: +31 (0) 50 528 04 98
 Fax: +31 (0) 50 528 04 95
 Mail: info@cobramedical.nl
 www.cobramedical.nl

Contact België

Tel: +32 (0) 70 66 05 70
 Fax: +32 (0) 70 66 05 72
 Mail: info@cobramedical.be
 www.cobramedical.be