


PEG

Used for per cutaneous endoscopic placement to provide enteral access for the delivery of nutrition, medication, gastric drainage or decompression



- Medical grade silicone with soft dome.
- Centimeter marking on the tube.
- Easy traction removal.
- Available in 20F & 24F with full kit version & economy version.
- Dual port feeding adapter available separately

Description	Order No.	QTY/BOX
PEG Kit 20 F	P-PEG 20	1
PEG Kit 24 F	P-PEG 24	1
PEG Kit 20 F Economy	P-PEG 20-E	1
PEG Kit 24 F Economy	P-PEG 24-E	1
Dual Port Feeding Adapter 20F	DFA 20	5
Dual Port Feeding Adapter 24F	DFA 24	5

Full Kit includes:

- 1 ea Silicone PEG feeding tube with soft internal dome
- 1 ea 120" insertion wire, 1 ea Scalpel with # 11 surgical blade
- 1 ea Dual Port Feeding Adapter, 1 ea 5cc disposable syringe
- 1 ea 25G x 1" needle, 1 ea 14 G needle/Cannula
- 1 ea External Bolster, 1 ea Drape, 1 ea Mosquito Forceps, 1 ea Clamp

Economy Kit includes:

- 1 ea Silicone PEG feeding tube with soft internal dome
- 1 ea 120" insertion wire
- 1 ea Scalpel with #11 surgical blade
- 1 ea Dual Port Feeding Adapter
- 1 ea External Bolster, 1 ea Clamp
- 1 ea 14g needle / cannula



For kit with retrieval snare add suffix - S