

1-Lumen Rotatable Sphincterotomes



RotaCut® – Sphincterotome with Braided Cutting Wire

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter
GSP-01-17-020	Standard Tip	5 mm	20 mm	1.7 mm
GSP-01-17-030	Standard Tip	5 mm	30 mm	1.7 mm
GSP-04-17-020	Tapered Tip	5 mm	20 mm	1.7 mm
GSP-04-17-030	Tapered Tip	5 mm	30 mm	1.7 mm



RotaCut® – Sphincterotome with Braided Cutting Wire

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter
GSP-02-17-020	Precut	-	20 mm	1.7 mm
GSP-02-17-030	Precut	-	30 mm	1.7 mm



Needle Knife Sphincterotomes

Catalog No.	Needle Projection	Outer Sheath Diameter
GSP-11-17-005	0 - 5 mm	1.7 mm
GSP-11-17-010	0 - 10 mm	1.7 mm
GSP-11-17-015	0 - 15 mm	1.7 mm

Special Versions



Special Sphincterotomes

Catalog No.	Description	Outer Sheath Diameter
GSP-12-18-004	Fistulotom, Cutting Spiral 4 mm	1.8 mm
GSP-13-17-015	Shark Fin Sphincterotome, Cutting Wire 15 mm	1.7 mm

Sphincterotomes

- Rotatable Tip  
*Precise Orientation & Positioning*
- 2 & 3 Lumen Versions  
*Contrast Medium Injection while Guide Wire is in Place*
- Monofilament Cutting Wire  
*Accurate & Accurate Cut*
- Version with Semi-insulated Cutting Wire  
*Safe & Controlled Papillary Roof Incision*
- Color Marking with Black Rings  
*Macroscopically Controlled Cannulation Depth*



Distributed by:

Ver. 4.0 (2010-08)

Manufactured by:  
Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany  
☎ +49 80 32/973-379 | 📠 +49 80 32/973-339 | ✉ sales@medi-globe.de

**Sphincterotomes**

Medi-Globe offers a wide range of sphincterotomes in triple, double and single lumen formats. A differentiating feature of Medi-Globe sphincterotomes is the ability to rotate the distal tip from 11 o'clock to 1 o'clock for precise orientation and positioning. Medi-Globe sphincterotomes are available with cutting wires of various lengths and in several tip styles to facilitate biliary and pancreatic duct access in nearly all anatomical requirements.

**RotaCut®**

RotaCut sphincterotomes are equipped with a braided cutting wire.

**RotaCut® - M**

RotaCut - M sphincterotomes are equipped with a monofilament cutting wire.


**Premier Cut**

Premier Cut sphincterotomes are equipped with a "semi-insulated", braided cutting wire.

**S.O.S. Sphincterotomes**

A S.O.S. "Speed Option Set" is comprised of a double lumen sphincterotome with braided cutting wire, which is preloaded with a hydrophilic Teflon coated Nitinol guide wire in .035" diameter and a length of 450cm.

**Sphincterotome Tips**



Medi-Globe sphincterotomes are directly compatible with Olympus electro-surgical generators. To verify compatibility with other electro-surgical generators please contact your Medi-Globe sales representative. Please also consult the Medi-Globe "Accessories" leaflet to see a complete list of active cords manufactured to fit most generators.

**3-Lumen Rotatable Guide Wire Sphincterotomes**



**RotaCut® – Sphincterotome with Braided Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Max. Accepted Guide Wire
GSP-31-20-020	Tapered Tip	5 mm	20 mm	2.55/2.05 mm	.035"
GSP-31-20-025	Tapered Tip	5 mm	25 mm	2.55/2.05 mm	.035"
GSP-31-20-030	Tapered Tip	5 mm	30 mm	2.55/2.05 mm	.035"

**RotaCut® - M – Sphincterotome with Monofilament Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Max. Accepted Guide Wire
GSP-32-20-020	Tapered Tip	5 mm	20 mm	2.55/2.05 mm	.035"
GSP-32-20-025	Tapered Tip	5 mm	25 mm	2.55/2.05 mm	.035"
GSP-32-20-030	Tapered Tip	5 mm	30 mm	2.55/2.05 mm	.035"

**Premier Cut – Sphincterotome with Braided and Semi-insulated Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Max. Accepted Guide Wire
GSP-27-20-030	Tapered Tip	5 mm	20 mm	2.55/2.05 mm	.035"

**2-Lumen Rotatable Guide Wire Sphincterotomes**



**RotaCut® – Sphincterotome with Braided Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Max. Accepted Guide Wire
GSP-21-20-020	Tapered Tip	5 mm	20 mm	2.05/2.55 mm	.035"
GSP-21-20-030	Tapered Tip	5 mm	30 mm	2.05/2.55 mm	.035"
GSP-24-20-020	Tapered Tip	10 mm (long nose)	20 mm	2.05/2.55 mm	.035"
GSP-24-20-030	Tapered Tip	10 mm (long nose)	30 mm	2.05/2.55 mm	.035"
GSP-23-20-030	Tapered Tip	20 mm (long nose)	30 mm	2.05/2.55 mm	.035"
GSP-21-17-020	Tapered Tip	5 mm	20 mm	1.7 mm	.021"
GSP-21-17-030	Tapered Tip	5 mm	30 mm	1.7 mm	.021"

**2-Lumen Rotatable Guide Wire Sphincterotomes**



**RotaCut® – Sphincterotome with Braided Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Max. Accepted Guide Wire
GSP-22-20-020	Precut Tip	-	20 mm	2.05/2.55 mm	.035"
GSP-22-20-030	Precut Tip	-	30 mm	2.05/2.55 mm	.035"
GSP-22-17-020	Precut Tip	-	20 mm	1.7 mm	.021"
GSP-22-17-030	Precut Tip	-	30 mm	1.7 mm	.021"

**RotaCut® – Billoth II Sphincterotome with Braided Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Max. Accepted Guide Wire
GSP-34-20-020	Billoth II	-	Variable	2.05/2.55 mm	.035"

**Premier Cut – Sphincterotome with Braided and Semi-insulated Cutting Wire**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Accepted Guide Wire
GSP-25-20-020	Tapered Tip	5 mm	20 mm	2.05/2.55 mm	.035"
GSP-25-20-030	Tapered Tip	5 mm	30 mm	2.05/2.55 mm	.035"
GSP-26-20-020	Tapered Tip	15 mm (long nose)	20 mm	2.05/2.55 mm	.035"
GSP-26-20-030	Tapered Tip	15 mm (long nose)	30 mm	2.05/2.55 mm	.035"

**S.O.S. Sphincterotome – "Speed Option Set" with Pre-loaded Hydrophilically Teflon Coated Nitinol Guide Wire in .035" Diameter and 450 cm Length**

Catalog No.	Description	Tip Length	Cutting Wire Length	Outer Sheath Diameter	Accepted Guide Wire
GSP-21-20-031	RotaCut®, Tapered Tip	5 mm	30 mm	2.55 / 2.05 mm	Pre-loaded
GSP-26-20-031	Premier Cut, Semi-insulated Cutting Wire, Tapered Tip	5 mm	30 mm	2.55 / 2.05 mm	Pre-loaded