

SonoTip® II EUS-FNA

Endoscopic Ultrasound-Guided FNA Systems – Single Use

SonoTip® II EUS-FNA


- Specially Treated Needle
Excellent Echogenicity
- Unique Sheath Stabilizer
Controlled Needle Advancement in Large Channel Scopes
- Lockable Sheath & Needle Length Adjustor
Variable Needle Projection (0-8.5 cm)
- Luer-locking Needle and Stylet
Secure Needles and Stylet Tip Positioning



SonoTip®II EUS-FNA


25 Gauge Needles

SonoTip®II – Including a Vacuum Syringe with a Filling Volume of 10 ml


Catalog No.	Description	Needle Diameter	Outer Diameter	Round Stylet
GUS-01-18-025	EUS-Guided FNA System → Plastic Outer Sheath	25 Gauge / 0.5 mm	1.8 mm	
GUS-01-27-025	EUS-Guided FNA System → Plastic Outer Sheath with Tip Stabilizer to Accommodate Large Working Channels (e.g. 3.2 mm, 3.7 mm, 3.8 mm)	25 Gauge / 0.5 mm	2.7 mm	

22 Gauge Needles

SonoTip®II – Including a Vacuum Syringe with a Filling Volume of 10 ml


Catalog No.	Description	Needle Diameter	Outer Diameter	Round Stylet
GUS-01-18-022	EUS-Guided FNA System → Metal Spiral → Compatible with Guide Wires of maximal .018"	22 Gauge / 0.7 mm	1.8 mm	
GUS-01-27-022	EUS-Guided FNA System → Metal Spiral with Tip Stabilizer to Accommodate Large Working Channels (e.g. 3.2 mm, 3.7 mm, 3.8 mm) → Compatible with Guide Wires of maximal .018"	22 Gauge / 0.7 mm	2.7 mm	

SonoTip®II – Including a Vacuum Syringe with a Filling Volume of 10 ml

Catalog No.	Description	Needle Diameter	Outer Diameter	Beveled Stylet
GUS-02-18-022	EUS-Guided FNA System → Metal Spiral → Compatible with Guide Wires of maximal .018"	22 Gauge / 0.7 mm	1.8 mm	
GUS-02-27-022	EUS-Guided FNA System → Metal Spiral with Tip Stabilizer to Accommodate Large Working Channels (e.g. 3.2 mm, 3.7 mm, 3.8 mm) → Compatible with Guide Wires of maximal .018"	22 Gauge / 0.7 mm	2.7 mm	


19 Gauge Needle

SonoTip®II – Including a Vacuum Syringe with a Filling Volume of 10 ml


Catalog No.	Description	Needle Diameter	Outer Diameter	Round Stylet
GUS-01-21-019	EUS-Guided FNA System → Metal Spiral → Compatible with Guide Wires of maximal .030" → (Recommended Medi-Globe catalog no. GGW-08-30-400)	19 Gauge / 1.0 mm	2.1 mm	

Accessories for SonoTip II EUS-FNA Line


Aspiration Syringe with a Larger Filling Volume of 20 ml

Catalog No.	Description	Aspiration Syringe
VAC120	VacLok Syringe 20 ml	

Nitinol Guide Wire Compatible to 19 Gauge EUS-FNA Systems


Catalog No.	Description	Diameter	Length	Nitinol Guide Wire
GGW-08-30-400	Nitinol Guide Wire, uncoated with 15 cm Radiopaque Flexible Spring Tip	.030"	400 cm	

Pseudo-Cyst Drainage Stents, Double Pigtail with Position Markers

Catalog No.	Description	Size	Length	Pseudo-Cyst Drainage Stent
GBS-11-08-002	Double Pigtail Stent with Position Markers	8.5 FR	2 cm	
GBS-11-08-005		8.5 FR	5 cm	
GBS-11-08-007		8.5 FR	7 cm	
GBS-11-08-009		8.5 FR	9 cm	
GBS-11-10-002		10.0 FR	2 cm	
GBS-11-10-005		10.0 FR	5 cm	
GBS-11-10-007		10.0 FR	7 cm	
GBS-11-10-009		10.0 FR	9 cm	

SonoTip[®] II EUS-FNA

Endoscopic Ultrasound-Guided FNA Systems – Single Use



Distributed by:

Ver. 1.0 (2010-10)

Manufactured by: